

Tooxy

A wide range of sensors

Voltmeter, ammeter, pH-meter, thermometer, light meter, hygrometer, pressure-meter, heart frequency, spirometry and ventilation, oximetry, CO2-meter



▲ Very simple heart frequency sensor with hand measurement.



▲ For physics, sensor like voltmeter or ammeter have two ranges for guaranteed results.



▲ Make 3 different measures with the ventilation sensor: flow, volume and respiratory frequency.

Tooxy

You plug in • You experiment



Designed by Jeulin,
exclusively engineered for science education



Rue Jacques-Monod, ZI n°1 Nétreville, BP 1900, 27019 Evreux cedex, France
Métropole • Tél : +33 (0)2 32 29 40 00 • Fax : +33 (0)2 32 29 43 99
International • Tél : +33 (0)2 32 29 40 23 • Fax : +33 (0)2 32 29 43 24
www.jeulin.com • export@jeulin.fr
SA au capital de 3 233 762 € - Siren RCS Evreux B 387 901 044 - Siret 387 901 044 00017

Article : 951288



91951288 10061010



www.jeulin.com

+33 (0)2 32 29 40 23



Measurement device or computer-assisted experimentation



Unique design



User friendly



Fish tape




USB port



Simplicity !

No need to install any driver or software:
Embedded in to the console, SCIENCE STUDIO runs automatically whatever your operating system.



Measurement device

Insert the sensor of your choice, the measured value appears on it.



Computer-assisted experimentation

Connect your Tooxy to a computer using a USB port, it is automatically detected and immediately operational.



Compatible operating systems

